



## AUDIOMETER 600 M



The audiometer **600 M** is particularly adapted to the healthcare professionals in the occupational, school and university medicine

# AUDIOMETER 600 M



**WARRANTY 3 YEARS (EXCEPT FOR THE HEADPHONES) – FRENCH DESIGN AND MANUFACTURING**

## TECHNICAL INFORMATION

■ **Type** : 4 – Screening and Control (according to EN 60645-1) Medical Marking EC 0459 - Class IIa  
■ **Output** : AC – Air Conduction – Light standard headphones  
■ **Frequencies range (11 values)** : 125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz ( +/- 1%)  
■ **Sound level** : from -10 HLdB to 100 HLdB  
■ **Increment of sound level** : 5 dB  
■ **Power supply** : isolated USB link ensuring the electrical safety (according to EN 60601-1)

■ **Stimulus** : continuous or pulsed sound  
■ **Minimum computer requirements** : 1 USB port, RAM : 1 Go, CPU : 1 GHz, Free space on hard disk : 50 Mo, Graphic resolution : 1024 x 600  
■ **Possible operating systems** : Windows 7, Windows Vista, Windows XP, Windows 8, Windows 10, Windows 11.  
■ **Dimensions audiometric box** : 92 x 150 x 28 mm  
■ **Packaging** : shock-proof briefcase 280 x 310 x 100 mm  
■ **Total weight** : 1.3 kg (160 g for the audiometric box only).

Electronica Medical offers a complete range of devices for screening the hearing losses, which are dedicated to the educational and occupational medicine.