



ELECTRONICA MEDICAL

by Electronica Technologies



AUDITEST



The manual audiometer **AudiTest** is achieving an **accurate screening** of hearing faculties. Convenient and performant, its size and ergonomics allow the Healthcare professional to use it with just one hand.

CE
0459

AUDITEST



Lightweight, portable, **manageable**. It can be used with just one hand.



Perfect **ergonomics** thanks to its quiet keys and its LCD display



- **Continuous or pulsed** sound
- **Frequencies** from 125 to 8000 Hz, output level range from -10 to 100 HLdB



Perfectly adapted to the **screening in schools**



- Device **designed and manufactured in France**
- **Warranty : 3 Years ***



End user manual : English, German, Italian, Spanish, Portuguese, Dutch, Hungarian, Polish, Romanian, French

WARRANTY 3 YEARS (EXCEPT FOR THE HEADPHONES) – FRENCH DESIGN AND MANUFACTURING

TECHNICAL INFORMATION

- **Type** : 4 – Screening and Control (according to EN 60645-1) Medical Marking EC 0459 – Class IIa
- **Output** : AC – Air Conduction
- **Frequencies range (11 values)** : 125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz (+/- 1%)
- **Sound level** : from -10 HLdB to 100 HLdB

- **Increment of sound level** : 5 dB
- **Power supply** : 2 battery sticks 1,5V -Type AA / LR06 – Supplied
- **Stimulus** : continuous or pulsed sound
- **Headphones** : light, comfortable
- **Dimensions of the audiometric box** : 79 x 122 x 33 mm
- **Packaging** : soft kit
- **Total weight of the kit** : 600 g (155 g for the audiometric box only).

Electronica Medical offers a complete range of devices for screening the hearing losses, which are dedicated to the educational and occupational medicine.